

AutoCPAP-20

CONTINUOUS POSITIVE
AIRWAY PRESSURE
(CPAP) THERAPY

Available for the rental basis and sale basis

CPAP devices are the most commonly used treatment options for people experiencing breathing problems, most notably of which is obstructive sleep apnea (OSA).

Product Specification

Characteristics	
Modes	CPAP, Auto-CPAP
Pressure range	4-20 CmH2O
Ramp time	0-60 minutes
Noise level	<28db
Alarms	Tube alarm, Apnea alarm & power failure
Memory	Micro SD Card - 8 GB
Heated Humidifier	Integrated
Warranty	1 Year

Physical Characteristics	
Dimensions	18*27*11 cm
Weight	1.7 Kg
Screen size	2.4 inch

Electrical Characteristics	
Electrical supply	100 - 240 V, 50 - 60 Hz
Power Consumption	≤ 80W with Humidifier
93/42/EC Classification	Class II A

Mode -> APAP			
Mode	APAP		
CPAP Start	4.0	cm H2O	
CPAP Max	11.6	cm H2O	
CPAP Min	4.0	cm H2O	
Ramp	5	min	
C-Flex	3		

Environmental requirements	
Temperature	5 C - 35°C
Humidity	30-95%, no condensation
Storage	-20-60-°C
Altitude	From Sea level to 2700 mtrs

Benefits of cpap

- Elimination of snoring and breathing obstructions.
- Improvement in quality of nightly sleep.
- Prevention or reversal of serious health conditions such as cardiovascular diseases and stroke.
- Lower blood pressure both during the day and at night.
- Increased alertness during the day.
- Significant decrease of daytime drowsiness.

2.1 Accessory List

No.	Name	Unit	Qty
1	APAP20/CPAP20 Device	Set	1
2	Humidifier(Optional)	Set	1
3	Air Tubing(Optional)	Set	1
4	Mask(Optional)	Set	1
5	DC Power Supply	Set	1
6	Micro SD Card(Optional)	Set	1
7	Canvas Bag	Set	1
8	Filter	Set	1
9	User Manual	Set	1



RESPIRCare

P/18