



BLOOD PRESSURE MONITOR

Specification :

Measurement Method	Oscillometric Method	
Measurement Range	Pressure	0mmHg ~ 300mmHg
	Pulse	30 to 180 Beats/Minute
Measurement Accuracy	Pressure	± 3mmHg
	Pulse	± 5%
Pressurization	Automatic Pressurization	
Memory	120 Memories in two Groups with Date and Time	
Function	Low Battery Detection	
	Automatic Power-Off	
Power Source	4 x AAA batteries	
Unit Weight	Approx. 326 g (11.5oz.) (Excluding Battery)	
Unit Dimensions	Approx. 139 × 88 × 43mm (5.5"x3.5"x1.7") (LxWxH)	

An ISO 9001 : 2008, CE 0197 Certified
With W H O Indicator



BLOOD PRESSURE MONITOR

Unit Dimensions	Approx. 139 × 88 × 43mm (5.5"x3.5"x1.7") (LxWxH)

An ISO 9001 : 2008, CE 0197 Certified
With W H O Indicator