



OXYGEN CONCENTRATOR

MODEL : MOC10L

Specification :

Technical parameter	
Model	MOC10L
Flow Rate	0-10 LPM
Purity	93% (±3%)
Outlet pressure	0.04-0.07 MPA (6-10PSI)
Noise level	<50dB
Power	AC220V , 50/60Hz (31Hz); AC 110V 50/60Hz (31Hz)
Power Consumption	< 880w
LCD Display	Starting times, operating pressure, present working time, accumulated working time, preset time from 10mins to 40hours
Alarm	Power failure alarm; high/low pressure alarm
Net Weight	27 Kgs
Size	365x375x600mm



OXYGEN CONCENTRATOR

MODEL : MOC5L

Specification :

Model	MOC5L
Flow Rate	0-5L/min
Purity	93 ± 3%
Outlet Pressure(Mpa)	0.04-0.07
Sound Level db	≤50db
Power	AC230V, 50Hz
Power Consumption	<300W
LCD Display	Operating pressure, Present Working Time, Accumulating Time, presetting Time from 10mins to 40hours
Alarm	Power failure alarm, High&Low pressure alarm Temperature alarm, Low purity alarm.
Net Weight	16 Kg
Size (mm)	365x375x600mm