

# PISTON COMPRESSOR PRO MINI NEBULIZER

Model No. SBS-PRO



# PISTON COMPRESSOR NEBULIZER

MCN - S600C



## Specification :

CERTIFIED: CE 0120

|                          |                                       |
|--------------------------|---------------------------------------|
| Power                    | 220 V 50 Hz                           |
| Power Consumption        | 180 VA                                |
| Medication Capacity      | 6 ml                                  |
| Particle Size            | Respirable Fraction<br>0.5 To 5um 80% |
| MMAD                     | 2M um                                 |
| Sound Level              | Below 55 dBA                          |
| Medication Bowl Capacity | 13ml( CC )                            |

|                             |                                                        |
|-----------------------------|--------------------------------------------------------|
| Average Nebulization Rate   | 0.2ml / min                                            |
| Compressor Pressure Range   | 30 to 36 Psi<br>( 210 to 250 KPa /<br>2.1 to 2.5 bar ) |
| Operating Pressure Range    | 8 to 16 Psi<br>( 50 to 100 KPa<br>0.5 to 1.0 bar )     |
| Liter Flow Range            | 6-10 lpm                                               |
| Operating Temperature Range | 50°F to 104°F<br>( 10°C to 40°C )                      |
| Operating Humidity Range    | 10 to 95 % RH                                          |

|                           |                                                                                   |
|---------------------------|-----------------------------------------------------------------------------------|
| Storage Temperature Range | -13°F to 158°F<br>(-21°C to 70°C )                                                |
| Storage Humidity Range    | 10 to 95 % RH                                                                     |
| Dimension ( L X W X H )   | 300 X 180 X 115 mm<br>( 11.8"X 7.1"X 4.5" )                                       |
| Weight                    | 4.5 Lbs ( 2.1kg )                                                                 |
| Standard Accessories      | Nebulizer, Air Tube ,<br>Mouth Piece, Filters<br>(5pcs) Adult Mask,<br>Child Mask |
| Optional Accessories      | Baby Pacifier Inhalator,<br>W/connector,<br>Nose Piece                            |