

LATEX FOLEY CATHETER - (SILICONISED COATED)



Specifications:

- Smooth tapered tip facilities for easy introduction thereby reducing urethral trauma.
- Symmetrical shaped balloon which expands equally in all directions.
- Accurately formed drainage lumen for effecient flow of urine.
- Larger circular drainage lumen for effecient flow of urine.
- Smooth flexible siliconised catheter shaft to provide patient comfort.
- Integrally bonded Reinforced junction between the inflation funnel and the catheter shaft to ensure reliability thus preventing detachment during use.
- Consistent shaft gauge.
- Plastic valves or Rubber valves fitted as per the convenience of the customers.
- Color circular sleeves for easy identification of size.
- Lateral opposing eyes (If requested by the customer).

Size	Balloon Capacity
2 Way 6 FR	3 ml - 5 ml
2 Way 8 FR, 10 FR	3 ml- 5 ml
2 Way 12 FR, 14 FR	15 ml - 30 ml
2 Way 16 FR, 18 FR, 20 FR, 22 FR, 24 FR	30 ml - 50 ml
3 Way 18 FR, 20 FR, 22 FR, 24 FR	30 ml - 50 ml

CE, ISO 9001 : 2008